

Kinematic Data Analysis using OpenCap

Presentation Date : 25.7.2


Presenter : Dongwoo Kim

What is OpenCap

OpenCap : Motion Capture analyzes Human Movement Dynamics from two smartphone videos




OpenCap : Low Cost & Less Time



Performance Validation of OpenCap

Knee loading during walking





Kinetic differences during sit-to-stand




Performance Validation of OpenCap

Muscle activation asymmetry during squatting




OpenCap Pipeline



OpenCap Experiment 1 : 2 Cameras


Motion : Squat



Simultaneous Video Capture



OpenPose



2D Keypoints


OpenCap Experiment 1 : 2 Cameras

Motion : Squat

LSTM


3D Keypoints

Triangulation



3D Anatomical Markers

OpenSim
IK Solver




Kinematics

OpenCap Experiment 1

Motion : Squat


Kinematics Data from the [OpenCap Website](#)



OpenCap Experiment 1


Motion : Squat

Kinematics Data from the OpenSim




OpenCap Experiment 1

Motion : Push-up




Simultaneous Video Capture



OpenCap Experiment 2 : 3 Cameras

Motion : Push-up

cam2




cam1




cam0



Simultaneous Video Capture






→ Accuracy has improved




Kinematics

Advantages of Using More than 2 Cameras



Advantages of Using More than 2 Cameras



Accuracy can be improved because Pose Estimation can be performed with various views.

Precautions when using OpenCap


It is advisable for the participant to wear **tight-fitting** clothing. Small Actions are not recommended.

Experiment 2 / squat / cam1



Joints are hard to detect.
Noise is generated by clothes.


Displacement of the **Left Knee**(3D Anatomical Markers)



Conclusion : Check the Data is Reliable

Check the OpenPose Results



Low Confidence Score Distribution




High Confidence Score Distribution




2D Keypoint & Confidence Score



Check the Filter Results before Dynamics



Displacement of 3D Anatomical Markers with Butterworth Filter



Thank You

Presenter : Dongwoo Kim